

Certificate

PROFIBUS Nutzerorganisation e.V. grants to

SIEMENS AG, DF FA AS DH AMB 3

Werner-von-Siemens-Str. 50, 92224 Amberg, Germany

the Certificate No: **Z02162** for the PROFIBUS Master:

Model Name: CPU1516TF-3 PN/DP
 Revision: 6ES7 516-3UN00-0AB0; SW/FW: V2.5; FS: 01
 GSD: si_m819c.gsd, File Version: 15.11.2017

This certificate confirms that the product has successfully passed the certification tests with the following scope:

<input checked="" type="checkbox"/> DP-V0	MS0, Sync, Freeze, Fail_Safe
<input checked="" type="checkbox"/> DP-V1	MS1, Prm_Block_Structure, Alarms
<input checked="" type="checkbox"/> DP-V2	IsoM
<input checked="" type="checkbox"/> Physical Layer	RS485

Test Report Number: M067-7
 Authorized Test Laboratory: SIEMENS AG, Fürth, Germany

The tests were executed in accordance with the following documents:
 "Test Specifications for PROFIBUS DP Masters, Version 3.0 from November 2005".

This certificate is granted according to the document:
 "Framework for testing and certification of PROFIBUS and PROFINET products".

For all products that are placed in circulation by **December 12, 2020** the certificate is valid for life.

Karlsruhe, June 18, 2018



(Official in Charge)

Board of PROFIBUS Nutzerorganisation e. V.



(Karsten Schneider)




(Dr. Jörg Hähnliche)